



SPECIFICATION AND DATA SHEET

RIGID PENDANTS



There are three basic types of MIPS Rigid Pendants:

MINI – available with up to 4 gasses and 4 Single Switched Sockets.

MIDI – available with up to 9 gasses and 8 Double Switched Sockets.

MAXI – as MIDI, but adjustable length by 300mm.

Each is designed for aesthetics, economy, maximum design life and reliability, whilst requiring minimum maintenance. They are available as a standard product with the following design features:

- ✓ Comply fully with the requirements of HTM 02-01
- ✓ Sturdy Monocoque design
- ✓ Manufactured from steel and epoxy powder coated
- ✓ Standard colour white
- ✓ All gasses fitted with BS 5682 N.I.S.T. hoses for easy installation
- ✓ First fix available separate and supplied with blanking facility
- ✓ Every Rigid Pendant pressure tested prior to despatch
- ✓ Every Rigid Pendant electrically tested prior to despatch

In addition to all of the above, MIPS can CUSTOM BUILD Rigid Pendants to your specific requirements, here are some examples:

- ✓ Other colours available
- ✓ Rigid copper pipework
- ✓ Special electrical connections (e.g. data, video, etc)
- ✓ Other makes of gas outlets available

Medical & Industrial Pipeline Systems Ltd, Dunshaughlin, Co. Meath, IRELAND

Tel: 00 353 1 825 0677 Fax: 00 353 1 825 0678 e-mail: sales@mips.ie web: <http://www.mips.ie>



MEMBER

