

Press Release August 2016

AREA VALVE SERVICE UNIT (A.V.S.U. Unit)

Medical and Industrial Pipeline Systems Ltd (MIPS) are proud to announce that they have commenced manufacture of their own Area Valve Service Units at their factory in Dunshaughlin, Co Meath.

Medical & Industrial Pipeline Systems Limited is the leading specialist in the design, manufacture, installation, commissioning & maintenance of piped gas equipment for the Medical, Industrial, Laboratory, Computer, Pharmaceutical, Educational and Electronic Industries in Ireland.

In-house manufacturing enables MIPS to the unit they can now tailor the A.V.S.U's to individual customer requirements.

The company offer a full one year warrants after installation and can provide an ongoing maintenance service.

Thus a customer can engage MIPS to design, install and commission full pressurised gases systems knowing that MIPS have been installing specialised gases equipment for over 40 years. MIPS are an ISO 9001:2008 company and member of British Compressed Gases Association.

SPECIFICATION AND DATA SHEET

The M.I.P.S. Area Valve Service Unit is designed for economy, maximum design life and reliability whilst requiring minimum maintenance. They are available as a standard product with the following design features:

CE marked

Powder coated steel enclosure

Fully removable quick release door for easy access

Comply fully with the requirements of HTM 02-01


All valves fitted with B.S.5682 N.I.S.T. connectors

Sealed shut off N.I.S.T. valves

All valves fitted with flat faced 'o' ring unions

Standard pipework sizes 15 – 28mm, others available upon request

Each valve supplied with blanking facility



Every A.V.S.U. assembly is pressure tested prior to despatch

Every valve tested for tightness (no pressure pass when closed)

In addition to all of the above, M.I.P.S. can offer the following options:

Specific RAL Colours

42mm pipework valves

Client specific units (size/depth/format)

For Further information contact



Medical & Industrial Pipeline Systems Ltd
Dunshaughlin Business Park
Dunshaughlin
Co. Meath
Ireland

Phone: 00353(0)18250677