


## LABORATORY GAS MANIFOLD RANGE

The MIPS range of Spectron manifolds are designed for economy, reliability and durability.



Pressure control panel BM56-2U



Pressure control panel BT2000AC-2x1-M-KH-SV-0 for acetylene with SV73



CRS2000V-2L automatic change-over model



CRS2000V-2 manual change-over model



CRS2000V-1 single-sided model

### Features and Benefits

- ✓ Fully complies with BCGA Codes of practice and covered by full warranty
- ✓ Robust wall mounted framework (floor mounted option available).
- ✓ Modular design leads to flexible extensions
- ✓ Highly visible gauges for cylinder contents and pipeline pressure.
- ✓ High pressure cylinder connection hoses to BSENISO14113.
- ✓ Cylinder support chains and racks included.
- ✓ Includes safety relief valve and isolation valve.
- ✓ Full range of models up to quality 6.0
- ✓ Each manifold is covered by a full warranty. Optional External weatherproof cages.

